



BD Influx Cell Sorter

Technical Specifications

The BD Influx™ cell sorter is a flexible flow cytometry platform that adapts to a researcher's application requirements. A modular architecture and a combination of detection capabilities, hands-on controls, and high performance allow researchers to configure the BD Influx system to various site and application needs.

The BD Influx system can handle a throughput rate of up to 200,000 events per second.

Optics

Excitation Optics

Optical Platform

Lasers are mounted on a standard optical breadboard. Laser beams are aligned independently by adjustable mirrors. Each laser beam has its own final focus lens mounted on a dedicated translational stage for fine adjustments. For six- and seven-laser systems, one or two laser paths have collinear beams through the steering optics and focusing lens, but intercept the stream independently.

Power Out of the Laser Head

355 nm: >100 mW
 405 nm: >50, >100 mW
 457 nm: >300 mW
 488 nm: >200 mW
 532 nm: >150 mW
 561 nm: >150 mW
 640 nm: >120 mW

Optional beam shaping optics allow for round or elliptical beam shapes (3:1 ratio with a typical beam height of 15–20 μm).

Emission Optics

Emission light is collected through a 20X, 0.6 NA microscope objective (90°). Light is focused on three, five, or seven spatially separated mirror pinholes depending on the number of lasers. Modular detector blocks allow for a user-defined detection configuration (see the *BD Influx Cell Sorter Filter Guide*).

Simultaneous video observation of the stream, the pinholes, and the laser intercepts allows for fast and intuitive alignment.

Regular forward scatter is detected using 75- and 50-mm lenses, an aperture, and a photomultiplier tube. Resolution for the standard forward scatter detector is >0.5 μm (measured using beads). Collection angle is 2–17°.

Side scatter is collected through the 90° collection lens and measured using a photomultiplier tube. Side scatter resolution is >0.2 μm (measured using beads and 0.1- μm filtered sheath fluid).

Fluidics

General Operation

- Laboratory air pressure and/or vacuum can be used for operation (regulated at 90 psi, 6.2 bar).
- An optional air pressure supply and/or vacuum pump is available.
- Sheath pressure is adjustable from 1–90 psi (0.07–6.2 bar).

Fluidics Reservoirs

Autoclavable 7-L sheath and waste containers, equipped with pressure and vacuum readout, are provided.

Fluidics Control

- Sheath, sample, and boost pressure can be individually adjusted.
- A sample flow fine adjustment is provided for precise regulation of sample flow.
- Purge, pulse, rinse, and run buttons are provided for quick stream startup and bubble removal.

Replaceable Fluidics Path

- The fluidics path, including the nozzle assembly, can be exchanged. There are no inline valves. Only pinch valves are used.
- The sample line can also be exchanged.

Bubble Detector

A bubble detector in the sample line detects air bubbles from the sample tube and stops sample flow when the sample tube is empty, preventing air bubbles from reaching the nozzle assembly.

Sample Input

12 x 75-mm tubes, polypropylene

Temperature Control

Sample tubes can be cooled or heated by an optional circulating water bath.

Performance

Fluorescence Sensitivity

Measured using SPHERO™ Rainbow Calibration Particles (RCP-30-5A) according to the manufacturer's specifications:

Sheath pressure: 33 psi
 Drop drive: On, ~68 kHz
 Excitation: 488 nm, 200 mW
 Emission: 530/40 nm for FITC
 580/30 nm for PE

FITC: 125 molecules of equivalent soluble fluorochrome (MESF-FITC)
 PE: 125 molecules of equivalent soluble fluorochrome (MESF-PE)

Fluorescence Resolution

Measured using propidium iodide (PI)-stained chicken erythrocyte nuclei (CEN):

Sheath pressure: 33 psi
 Drop drive: On, ~68 kHz
 Excitation: 488 nm, 200 mW
 Emission: 610/20 nm for PI

Coefficient of variation (CV) of PI: <3%, full G_0/G_1 peak

Fluorescence Linearity

Measured using propidium iodide PI-stained CEN:

Sheath pressure: 33 psi
 Drop drive: On, ~68 kHz
 Excitation: 488 nm, 200 mW
 Emission: 610/20 nm for PI

Doublet/singlet ratio: 1.95–2.05

Sort Performance

Drop Drive Frequency

Adjustable 9–180 kHz

Purity and Yield

At 60 psi and 100 kHz with an average threshold rate of 25,000 events per second, a four-way sort achieved a purity of 98% and a yield >80% of Poisson's expected yield for all four populations. Higher threshold rates of up to 70,000 events per second can be achieved without affecting purity; however, yield will decrease based on Poisson statistics.

Viability

As shown in published literature, sorts performed using murine¹⁻⁴ and human⁵ cells and/or cell lines⁶ demonstrated good recovery and viability in several experimental systems. Optimal sort conditions need to be established for different cell types.

Nozzles

Supplied: 70, 86, 100, and 140 μm

Optional: 200 μm

Sort Collection Devices

All collection devices are designed to fit on the Computerized Cell Deposition Unit (CCDU). The CCDU is standard on all instruments.

- Two-way sorting: Microtubes; 12 x 75-mm, 15-mL, and 50-mL tubes
- Three-way sorting: One 50-mL tube and two 12 x 75-mm tubes
- Four-way sorting: Microtubes, 12 x 75-mm tubes
- Six-way sorting:* 0.4-mL tubes
- Plates and slides: 6, 24, 48, 96, and 384-well plates; slides; and user-defined collection devices

Sample Collection Cooling

Water recirculator for refrigeration or heating of collection tubes (optional).

Sort Monitoring

- Live video feed of waste collection and side streams.
- Live video feed of breakoff point.
- Drop-delay determination is achieved with BD FACS™ Accudrop technology. The drop-delay value can be adjusted while viewing BD Accudrop beads which are illuminated by a red diode laser in the center and side sort streams.

Signal Processing

Data Acquisition Channels

- Up to 5-laser systems: 16 channels, usually 14 colors plus forward and side scatter
- 6- and 7-laser systems: 24 channels,* usually 22 colors plus forward and side scatter

Signal Processing

- 16-bit analog-to-digital conversion, 65,536 channels
- Parallel data stream with channel ID and integrity check
- Less than 1 correlation error per 10^8 events

Acquisition Rate

Dead time is 0 μs . The maximum throughput rate is 200,000 events per second, independent of the number of parameters.

Fluorescence Compensation

- Up to 5-laser systems: 16 x 16 digital compensation matrix.
- 6- and 7-laser systems: 24 x 24 digital compensation matrix.*
- Compensated parameters are being added to the bus as separate parameters.

Pulse Processing

- All signals are height (peak) by default.
- Optional pulse processor electronics add area and width measurements for a maximum of 8 parameters to the bus.
- Width measurement on the trigger parameter is standard.

Time

Time can be correlated to any parameter for kinetic experiments or other applications.

Channel Threshold

- Any parameter can be used as the threshold from the primary laser.
- Lasers and detectors can easily be switched to change laser sequence.

Workstation

Workstation

PC workstation with at minimum Intel® Core2 Duo, 2.40 GHz or faster, USB (2 ports), Windows® XP professional operating system

Memory

1 GB

Data Storage

- 2 x 80 GB hard drives, RAID 1 (mirrored) configured
- 16x CD/DVD +/- RW
- Floppy drive

Networking

10/100/1000 Ethernet

Monitor

20-inch LCD, 1280x1024 resolution

Data File Structure

Flow cytometry standard (FCS) 3.0

Software

BD Spigot v6.1.4 or later

* Enabled with new software currently in development.

Installation Requirements

Dimensions (H x D x W)

BD Influx: 55.1 x 50.4 x 39.4 in.
(140 x 128 x 100 cm)

BD Influx with HEPA-filtered enclosure:*
85.8 x 50.4 x 49.2 in.
(218 x 128 x 125 cm)

BD Influx with HEPA-filtered enclosure
(extended bench):*
85.8 x 50.4 x 59.1 in.
(218 x 128 x 150 cm)

Weight does not exceed 485 lb (220 kg).

* A minimum of 102.4 in. (260 cm) floor-to-ceiling height is required during installation of the HEPA-filtered enclosure.

Size and Weight

On all sides, 18 in. (46 cm) is required for service access. A minimum of 102.4 in. (260 cm) floor-to-ceiling is required for installation with the HEPA-filtered enclosure.

Temperature Operating Range

18–22 ±3°C

Heat Dissipation

6,200 BTU/h maximum, dependent on laser choice

Power

Two (2) dedicated 120-VAC, 50/60-Hz, 15-A lines. Multiple peripherals may require an additional line. An appropriate transformer will be delivered depending on region-specific settings. Optional air and vacuum supply require an additional line. Power requirements are between 1,080 watts per hour (3 lasers, no HEPA-filtered enclosure) and 1,800 (5 lasers with HEPA-filtered enclosure). Standby status is not applicable.

Humidity

55% ±10% RH

Room Noise

<80 dBA from all running equipment

Water Supply

Not required unless water-cooled lasers are used. Specifications defined by laser manufacturer.

Air Supply

90 psi (6.2 bar) regulated. The source of compressed air must deliver clean (less than 5 ppm), dry-filtered (oil-free) air and stable pressures.

Vacuum Supply

A vacuum supply is delivered with the system. A laboratory vacuum supply between 5 and 15 in. Hg at 1 CFM can also be used.

Options

Small particle

Small particle with polarization

HEPA-filtered enclosure

Aerosol evacuation

Sample temperature control

Polarization (without small particle)

Air compressor

Regulatory Status

CE marked for electrical safety (Europe)

UL Standard for electrical safety (USA)

CSA for electrical safety (Canada)

Class I (1) laser product per CDRH regulations and EN/IEC 60825-1

References

1. Chen HL, Pistollato F, Hoepfner DJ, Ni HT, McKay RDG, Panchision DM. Oxygen tension regulates survival and fate of mouse central nervous system precursors at multiple levels. *Stem Cells*. 2007;25:2291-2301.
2. Kumar R, Fossati V, Israel M, Snoeck HW. Lin-Sca1+kit-bone marrow cells contain early lymphoid-committed precursors that are distinct from common lymphoid progenitors. *J Immunol*. 2008;181:7507-7513.
3. Mar L, Hoodless PA. Embryonic fibroblasts from mice lacking Tgfr were defective in cell cycling. *Mol Cell Biol*. 2006;26:4302-4310.
4. Seamons A, Percheret A, Gorman J. Endogenous myelin basic protein is presented in the periphery by both dendritic cells and resting B cells with different functional consequences. *J Immunol*. 2006;177:2097-2106.
5. Pistollato F, Chen HL, Schwartz PH, Basso G, Panchision DM. Oxygen tension controls the expansion of human CNS precursors and the generation of astrocytes and oligodendrocytes. *Mol Cell Neurosci*. 2007;35:424-435.
6. Panchision DM, Chen HL, Pistollato F, Papini D, Ni HT, Hawley TS. Optimized flow cytometric analysis of central nervous system tissue reveals novel functional relationships among cells expressing CD133, CD15, and CD24. *Stem Cells*. 2007;25:1560-1570.

Class I (1) laser product.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Windows is a registered trademark of Microsoft Corporation.

Intel is a registered trademark of Intel Corporation.

SPHERO is a trademark of Spherotech, Inc.

BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. © 2010 BD



BD Biosciences
2350 Qume Drive
San Jose, CA 95131
answers@bd.com
bdbiosciences.com